

OIML Standard Weights – Class E1



E1

SPECIFICATION:

Class	E1
Standard	OIML R111
Material	Austenitic Stainless Steel
Construction	Solid single piece fabrication
Cavity	No cavity adjustment
Density	8.00g/cm³
Magnetic Susceptibility	< 0.005
Intensity of polarization	< 2.5 μ T
Shape	Cylindrical with knob
Finish	Mirror Polish
Packing	Aluminum alloy box
Accessories	Tweezer, Cloth, Gloves

Our Range Class E1

E1 - 1mg to 500mg - 12pcs. E1 - 1mg to 200g - 23pcs.

Application

E1 Class weights are recommended for primary reference standard for calibration other
It cannot use in General Laboratory.



ELS Instruments

**Call Us****+91 94488 27690****Reach Us at****enquiry@elsinstruments.com**