


Chemical Storage



CHEMSTOR-CHEMICAL STORAGE CABINETS

Since hazardous materials must be stored safely, Els Instruments manufactures chemical storage cabinets with no exposed metal parts. Units are constructed out of acid resistant polypropylene. All joints are thermally welded by a special process without using metal fasteners.

Chemical storage cabinets can be connected to an extraction system for continuous extraction of volatile gases. This is to make sure that the laboratory air is safe to breathe and environment friendly.



ELS Instruments



Call Us
+91 94488 27690



Reach Us at
enquiry@elsinstruments.com