

## Sink and Drip Cups



Lab Sinks and Drip cups shall be selected to match the worktop and application. The varieties of sinks that Els Instruments offers are – Stainless steel sink, Polypropylene sink, Epoxy sink and Ceramic sink.

### **Description:**

Sinks	Suitable worktops
Polypropylene sink	Polypropylene, Epoxy, Phenolic resin, Acrylic Melamine laminate, Compact laminate
Epoxy sink	Epoxy, Phenolic resin
Stainless steel sink	Stainless steel, Melamine laminate
Ceramic sink	Ceramic worktops