

## Laboratory Modular Furniture



we Manufacturer and Supplier **Laboratory Modular Furniture** that is appreciated for reasonable pricing and innovative designing. Our research and development team takes individual requirements of the customers into account to provide them with user-friendly Laboratory Modular Furniture having a matchless perfection.



ELS Instruments has the provision of a wide variety of stylish lab furniture, providing modern work benches. We have different types of work benches varying in colour and designs – Wall bench, Island bench and Peninsular bench – each of which has a specific significance and are fit for a variety of purposes.

## ELS Basic Series

The worktops are installed on under bench cabinets–panels. The cabinets are constructed on integrated plinths of 100mm to keep the cabinet bottom elevated from the floor. The cabinet rests on the floor without any additional support of framework

## ELS Suspended Series

Worktops are mounted on heavy duty cantilever metal frame. Suspended cabinets or mobile cabinets shall be used with this system. Thus, under-bench spaces can be organized by the user to adapt their needs easily without the help of a technician. Ease of access to the service area is a special advantage of this series

## ELS H Series

Worktops are placed on metal frameworks. This series is designed to support heavier loads. The 'H' frames are made from mild steel square sections. Under bench cabinets rests on frames.

