


## **Work Tops**



Worktops are the most used part of a laboratory furniture. Hence, it shall be selected carefully to suit the application, chemical resistance, thermal stability, impact resistance, durability, ease of sterilization, aesthetics and cost effectiveness.

### **We have various types of worktops available such as:**

- Phenolic Resin ( Trespa )
- Epoxy Resin ( Durcon )
- Melamine laminated plywood
- Compact laminate
- Polypropylene
- Stainless steel
- Acrylic
- Ceramic worktops
- Granite
- Solid timber