

Analytical Balance



Analytical Balance – Touch Screen




MAB 220T

- Electro Magnetic Force Compensation
- Touch Screen display
- Percentage Weighing
- Multiple Weighing Units
- Density Determination Function
- RS232C Interface
- Counting Weighing mode
- Underhook Weighing
- Calibration with Internal weight
- Time
- Date

PRODUCT SPECIFICATIONS:

MODEL	MAB 220T (Touch screen)
Max. (Capacity)	220 g
Readability (d)	0.0001 g
Repeatability (+)	0.2 mg
Linearity (+)	0.4 mg
Min. (Capacity)	10 mg
Verification Scale Intervals (e)	1 mg
Calibration	Motorized Internal Calibration
Class	Class I
Display	5 inch Touch Screen
Response Time	4 - 5 Sec
Pan Size (mm)	80 mm
Tare Range	Full
Operating Temperature	+15 C to + 25 C
Sensitivity Drift +	2ppm / C
Density	Function Available
Key Sound	Yes

Weighing Chamber	240 x 190 x 265 mm
Unit	g, mg, ct, oz
Baud rate Options	1200, 2400, 4800, 9600
Interface	RS232C
Power Supply	AC 220V / 50-60 Hz
Dimension (LxWxH) mm	475 x 315 x 485
Weight	10 Kg



ELS Instruments

**Call Us****+91 94488 27690****Reach Us at****enquiry@elsinstruments.com**