

Autoclave Horizontal cylindrical



Triple walled unit made of S.S. 304 grade all joints are argon welded. Mounted horizontal on a sturdy. heavy mild steel tubular stand duty powder coated. Pressure is adjustable horn 10-20 PSI. Fitted with pressure gauge. Safety valve and steam release valve. The units fitted with MULTIPLE OPERATING VALVE to carry out all functions of sterilization cycle at one point i.e. Steam to chanter.

Exhaust, Vacuum/Drying. Fitted with a Digital Temperature Indicator to indicate inner temperature. Automatic Pressure Control device and Automatic Low Water Level Cut-off Device. Lid is made of S.S. plate provided with radial lacking system.

Triple walled unit made of S.S. 304 grade all joints are argon welded. Mounted horizontal on a sturdy. heavy mild steel tubular stand duty powder coated. Pressure is adjustable horn 10-20 PSI. Fitted with pressure gauge. Safety valve and steam release valve. The units fitted with MULTIPLE OPERATING VALVE to carry out all functions of sterilization cycle at one point i.e. Steam to chanter.

Exhaust, Vacuum/Drying. Fitted with a Digital Temperature Indicator to indicate inner temperature. Automatic Pressure Control device and Automatic Low Water Level cut-off Device. Lid is made of S.S. plate provided with radial lacking system.

Model No.	Dia x Depth	Capacity	Load
EAUC-HC1	250 x 450MM	22Ltrs	2.0KW
EAUC-HC2	300 x 500mm	35 Ltrs	3.0KW
EAUC-HC3	350 x 550mm	50 Ltrs	4.0KW
EAUC-HC4	400 x 600mm	75 Ltrs	6.0KW
EAUC-HC5	450 x600mm	95 Ltrs	6.0KW
EAUC-HC6	550 x750mm	175 Ltrs	8.0KW
EAUC-HC7	500 x1100mm	215 Ltrs	9.0KW

NB: Special sizes can be fabricated on request

Hot plate Rectanguar With Energy Regulator

Heavy MS Body duty powder coated with SS or MS Top, fitted oath Energy Regulator and Rotary Switch supplied complete with cord and plug to work on 220/230V AC (Temperature 250°C-300°Capprox.)

S.S Top

Model No.	Size	Watts
EHPE-1SS	250mmx300mm	2000
EHPE-2SS	250mm x400mm	2000
EHPE-3SS	300mm x450mm	2000
EHPE-4SS	450mm x600mm	3000

M.S Top

Model No.	Size	Watts
EHPE-1MS	250mmx300mm	2000
EHPE-2MS	250mm x400mm	2000
EHPE-3MS	300mm x450mm	2000
EHPE-4MS	450mm x600mm	3000