

# **Autoclave Triple Walled Vertical**



## **With Radial Locking System Autoclave**

Dry type vertical triple walled construction. all the walls are made of thick high grade stainless Vertical steel sheet of SS-304 grade. The lid is made of Stainless-steel lid with radial locking system and foot paddle lifting device with silicon gasket. provided with built-in safety valve, pressure gauge. pressure release valve and water level indicator. The pressure is adjustable from 5 psi with automatic pressure control Switch. controlled with microprocessor based digital controller with timer supplied with stainless steel basket.

Model No.	Dia x Height	Rating	Volume
EAUT-S1	250 x 450mm	1.5KW	22Ltr
EAUT-S2	300 x 500mm	2.0KW	35Ltr
EAUT-S3	350 x 550mm	3.0KW	53Ltr
EAUT-S4	450 x 600mm	4.0KW	90Ltr
EAUT-S5	550 x 750mm	6.0KW	175Ltr

## **OPTIONAL**

Lid with Bakelite knobs can be (instead of radial locking system) On request

- Auto Pressure cut off system
- Auto Vacuum Breaker
- Auto Purging & Exhaust (Solenoid Value)

**NB:** Special size can be fabricated on request