


(Borosilicate Boiler & Borosilicate Condenser & Quartz Heater with built in Triple Safety Feature).



TECHNICAL SPECIFICATIONS:

CAT. NO:	EBS- 2	EBS- 4	EBS- 8
Output	2 litres/hr, single distilled	4 litres/hr, single distilled	8 litres/hr, single distilled
Boiler	Borosilicate	Borosilicate	Borosilicate
Condensor	Borosilicate	Borosilicate	Borosilicate
pH	5.5 - 6.0	5.5 - 6.0	5.5 - 6.0
Conductivity, μScm^{-1}	1.0-2.0	1.0-2.0	1.0-2.0
Temperature	30 - 40°C	30 - 40°C	30 - 40°C
Pyrogen content *	Pyrogen free	Pyrogen free	Pyrogen free
Water supply	2 litre/min,	2 litre/min,	4 litre/min,
Electricity supply	220 or 240V, 50 -60Hz, Single phase	220 or 240V, 50 -60Hz, single phase	220 or 240V, 50 -60Hz, double phase
Power requirement	1.5kW	3kW	6kW

The above output result can be achieved by using treated water as a boiler feed.