

# **(Borosilicate Single Distillation) with Metal Heater**



## **TECHNICAL SPECIFICATIONS:**

<b>Output</b>	<b>4 litres/hr, single distilled</b>
pH	5.0-6.5
Conductivity, $\mu\text{Scm}^{-1}$	3.0-4.0
Resistivity, $\text{mOhm-cm}$	0.25-0.3
Temperature	25-35°C
Pyrogen content*	Pyrogen free
Water supply	1 litre/min,
Electricity supply	220 or 240V, 50-60Hz, single phase
Power requirement	3kW
Dimensions, (w x d x h)	500 x 150 x 450mm

The above output result can be achieved by using treated water as a boiler feed.