

(Borosilicate Single Distillation Basic model) with Quartz Heater



TECHNICAL SPECIFICATIONS:

CAT.NO:	EBSD-1.5	EBSD-4	EBSD-8
Output	1.5 litres/hr, single distilled	4 litres/hr, single distilled	8 litres/hr, single distilled
pH	5.0 - 6.5	5.0 - 6.5	5.0 - 6.5
Conductivity, μScm^{-1}	1.5-2.0	1.5-2.0	1.5-2.0
Temperature	30 - 40°C	30 - 40°C	30 - 40°C
Pyrogen content	Pyrogen free	Pyrogen free	Pyrogen free
Water supply	2 litre/min, 3 - 100psi, (20-700kPa)	3 litre/min, 3 - 100psi, (20-700kPa)	6 litre/min, 3 - 100psi, (20-700kPa)
Electricity supply	220or240V,50-60Hz, single phase	220or240V,50-60Hz, single phase	220or240V,50-60Hz, double phase
Power requirement	1.5kW	3kW	6kW

The above output result can be achieved by using treated water as a boiler feed.