

Deep Freezer Horizontal (Chest Type)



This is an ecofriendly and CFC free refrigerant and constructed of double wall, the exterior is made of MS Sheet. duly powder coated. interior of SS having PUF insulation to minimize heat loss. Capillary controlled, balanced and hermetically sealed refrigeration is most economical, trouble free and exceptionally quiet performance. Temperature is controlled by Micro Processor based Digital Temperature Controller with single compartment having an accuracy of $\pm 3^{\circ}\text{C}$

Upto -20°C

Model No.	Inn. Dimension [W x D x H]	Max Temp	Vol
EDF-1	500mmx380mmx600 mm	-20°C	115Ltrs
EDF-2	700mmx480mmx600 mm	-20°C	200Ltrs
EDF-3	1000mmx500mmx600 mm	-20°C	300Ltrs
EDF-4	1000mmx600mmx700 mm	-20°C	420Ltrs

Upto -35°C

EDF-1SP	500mmx380mmx600 mm	-35°C	115Ltrs
EDF-2SP	800mmx380mmx600 mm	-35°C	180Ltrs
EDF-3SP	1200mmx400mmx600 mm	-35°C	300Ltrs

OPTIONAL Automatic Voltage Stabilizer 3 Kva

NB: Non use of stabilizer will render the warranty of the equipment invalid