

Deep Freezer Vertical (Upright)



This is an ecofriendly and CFC free refrigerant and contacted of double wall, the exterior is made of M. S Sheet, duly powder coated and interior of S.S having PUF insulation to minimize heat loss. Capillary controlled. balanced hermetically sealed refrigeration is most economical, trouble free and exceptionally quiet performance and supplied with four Nos. standard storing trays. Temperature is controlled by Micro Processor based Digital Temperature Controller Cum Indicator having an accuracy of $\pm 3^{\circ}\text{C}$.

Upto-20°C

Model No.	Inn. Dimension [W x D x H]mm	Max Temp	Vol
EDFV-1	400x425x600	-20°C	115Ltrs
EDFV-2	500x425x800	-20°C	171Ltrs
EDFV-3	500x580x980	-20°C	285Ltrs
EDFV-4	650x640x1000	-20°C	425Ltrs

Upto-35°C

EDFVSP-1	400x425x600	-35°C	115Ltrs
EDFVSP-2	500x425x800	-35°C	171Ltrs
EDFVSP-3	500x580x980	-35°C	285Ltrs

OPTIONAL:Automatic Voltage Stabilizer 3 Kva

NB:Non use of stabilizer will render the warranty of the equipment invalid