

## **Flash Autoclave (Dental)**



It is tabletop front loading autoclave. double walled unit. outer M.S. powder coated inner 304 grade S.S. All joints are argon welded. Door is made up of heavy gauge S.S. sheet and using high quality silicon gasket used with ISI marked immersion heater 1.5 KW In 230/240 V.

### **Features:**

- Fully automatic microprocessor-based control
- Digital timer both for wet and dry cycle
- Alarm after the total command is completed
- Cyclic timer 10 to 15 minutes
- Temp range 121°C-134°C

Unit is supplied with 2 Nos. S.S. trays with tray stand

Model No.	Dia x Depth	Volume
EFAC-1	230 x 410mm	17Ltr

### **Autoclave Portable**

**Single** wall electrically operated constructed with Stainless steel 304 grade

Model N	Dia x Height	Body	KW	Capacity
EAUC-PS	300 x300mm	S.S	1.5	20Ltr