

Platelet Incubator cum Agitator



These platelet Incubators are specially designed for close environment control for stored platelets. Ideally suitable for accommodating ELS Platelet Agitator. Unit also consist of Platelet Agitator and working in a fixed temperature at 22°C having an accuracy of $\pm 0.5^\circ\text{C}$. Constructed of double walls. the exterior is made of mild steel duly powder coated, and the inner chamber is made of SS. Theplex glass inner door make it possible to inspect without disturbing the temperature. All controls are mounted at the bottom of Platelet Incubator.

Model No	Inner size	Capacity
EPI-110	480 x 380 x 580 mm	110Ltr

OPTIONAL

- a) Recorder 7 days circular chart Recorder- Inkless and ink type
- b) Automatic Voltage Stabilizer 1 KVA

NB: Nonuse of stabilizer will render the warranty of the equipment invalid