

with Glass Windshield
Precision Gold Balance
PGB 610



Salient Features

- Glass Windshield
- Auto Calibration with External Weight
- Low Battery Indicator
- Multifunction Weighing Units
- In built rechargeable battery for continuous use
- Auto Zero tracking
- Overload Protection Design
- Front RED LED display
- Automatic Switch Off Option
- High precision loadcell sensor
- RS232 Interface

SPECIFICATIONS

MODEL	PGB 630	PGB 1030
Capacity	600 g	1000 g
Readability	0.01 g (10 mg)	
Repeatability (±)	0.01 g	
Linearity (±)	0.02 g	
Pan Size	120 Ø mm	
Calibration	Automatic External	
Tare Range	Full	
Operating Temperature	10oC to 40 oC	
Power Supply	AC 220V / 50-60 Hz	
Dimension (LxWxH) mm	325 x 260 x 350	
Weight	3.7 kg	