

Precision Balance

PGB 3010 & 6000




Salient Features

- Auto Calibration with External Weight
- Easy to read LCD Display with backlight
- Multifunction Weighing Units
- Auto Zero tracking, Low Battery Indicator
- Overload Protection Design
- S.S. Weighing Pan, Level indicator, Adjustable Fit
- High Resolution, Quick Weighing, Accurate Result
- RS232C Interface
- In built rechargeable battery for continues use.

MODEL	PGB 630	PGB 1030
Capacity	3000 g	6000 g
Readability	0.01 g (10 mg)	0.1 g (100 mg)
Repeatability (±)	0.01 g	0.1 g
Linearity (±)	0.02 g	0.2 g
Pan Size	180 x 195 mm	

Display	LCD
Calibration	Automatic External
Tare Range	Full
Operating Temperature	10oC to 40oC
Power supply	AC 220V/ 50-60 Hz
Dimension (LxWxH) mm	440 x 265 x 220
Weight	4 kg



ELS Instruments

**Call Us****+91 94488 27690****Reach Us at****enquiry@elsinstruments.com**