

Precision Balance

PGB 3110 & 6100



Salient Features

- Auto Calibration with External Weight
- Easy to read LCD Display with backlight
- Multifunction Weighing Units
- Auto Zero tracking, Low Battery Indicator
- Overload Protection Design
- S.S. Weighing Pan, Level indicator, Adjustable Fit
- High Resolution, Quick Weighing, Accurate Result
- In built rechargeable battery for continues use.

SPECIFICATIONS

MODEL	PGB 630	PGB 1030
Capacity	3000 g	6000 g
Readability	0.01 g (10 mg)	0.1 g (100 mg)
Repeatability (±)	0.01 g	0.1 g
Linearity (±)	0.02 g	0.2 g
Pan Size	180 x 195 mm	
Display	LCD	
Calibration	Automatic External	
Tare Range	Full	
Operating Temperature	10oC to 40oC	
Power supply	AC 220V/ 50-60 Hz	
Dimension (LxWxH) mm	440 x 265 x 220	
Weight	4 kg	