

Ultra-Low Freezers Horizontal (chest type)



Up To -80 0 C [in Ac Room]

This freezer specially engineered superior insulation and specially designed refrigeration system- double walled with exterior of MS duty powder coated interior of SS having PUF insulation to minimize heat loss. Capillary controlled. hermetically sealed CPC free refrigeration system to ensure most economical / trouble free and exceptionally quiet performance and provided with castor wheels. Temperature is controlled by PID Micro Processor Based Digital Temperature Controller Cum Indicator [LED] with single compartment upto -80 -°C in AC room having an accuracy of +/-3 °C.

Model No.	Inn. Dimension [W x D x H]	Vol
EDFH-U1	550x370x680mm	110Ltrs
EDFH-U2	800x370x800 mm	170Ltrs
EDFH-U3	1200x370x980 mm	285Ltrs

Same as above but Temperature -40°C To-70°C [in AC room]

EDFH-L1	550x370x680mm	110Ltrs
EDFH-L2	800x370x800 mm	170Ltrs
EDFH-L3	1200x370x980 mm	285Ltrs

OPTIONAL Automatic Voltage Stabilizer 3 Kva

NB: Non use of stabilizer will render the warranty of the equipment invalid