

Ultra lowfreezers vertical(Upright)



Upto -80 °C [in AC room]

This freezer specially engineered superior, insulation and specially designed refrigeration system- double walled with exterior of MS with duly powder coated interior of SS having PUF insulation to minimize heat lost Capillary controlled. hermetically sealed CFC free refrigeration system to ensure most economical trouble free and exceptionally quiet performance. The unit is supplied with four no's standard storing trays with castor wheels. Temperature is controlled by PID Micro Processor Based Digital Temperature Controller [LED] up to -80 -C in AC room having an accuracy of + / -3°C.

Model No.	Inn. Dimension [W x D x H]	Vol
EDFV-U1	425x400x680mm	110Ltrs
EDFV-U2	500x425x845mm	170Ltrs
EDFV-U3	580x500x990 mm	285Ltrs

Same as above but Temperature -40°C To-70°C [in AC room]

EDFV-U1	425x400x680mm	110Ltrs
EDFV-U2	500x425x845mm	170Ltrs
EDFV-U3	580x500x990 mm	285Ltrs