

GLASS DOUBLE WATER STILLS **(Borosilicate Double Distillation Basic model) with Quartz Heater**



TECHNICAL SPECIFICATIONS:

CAT.NO:	EPTI-D-1.5	EPTI-D-4	EPTI-D-8
Output	1.5 litres/hr, Double distilled	4 litres/ hr, Double distilled	8 litres/hr, Double distilled
Boiler	Borosilicate	Borosilicate	Borosilicate
Condensor	Borosilicate	Borosilicate	Borosilicate
pH	5.0 - 6.5	5.0 - 6.5	5.0 - 6.5
Conductivity, μScm^{-1}	0.8-1.2	0.8-1.2	0.8-1.2
Temperature	30 - 40°C	30 - 40°C	30 - 40°C
Pyrogen content *	Pyrogen free	Pyrogen free	Pyrogen free
Water supply	60 litre/hr	90 litre/hr	120 litre/hr
Powerrequirement	3Kw	6kW	12kW

The above output result can be achieved by using treated water as a boilerfeed.