

(Double, Borosilicate Boiler & Borosilicate Condenser & Quartz Heater with built in triple safety feature).



TECHNICAL SPECIFICATIONS:

CAT.NO:	EBB-2	EBB- 4	EBB-8
Output	2 litres/hr, Double distilled	4 litres/hr, Double distilled	8 litres/hr, Double distilled
Boiler	Borosilicate	Borosilicate	Borosilicate
Condensor	Borosilicate	Borosilicate	Borosilicate
pH	6.0-7.0	6.0-7.0	
Conductivity, $\mu\text{Scm-l}$	0.8-1.0	0.8-1.0	0.8-1.0
Temperature	30 - 40°C	30 - 40°C	30 - 40°C
Pyrogen content *	Pyrogen free	Pyrogen free	Pyrogen free
Water supply	1.5 litre/min x 2	2 litre/min x 2 Input	4 litre/min x 2 Input
Electricity supply	220 or 240V,50-60Hz, single phase	220 or 240V,50-60Hz, single phase	220 or 240V,50-60Hz, single phase
Powerrequirement	3kW	6kW	12kW

The above output result can be achieved by using treated water as a boiler feed.