

(Double Quartz Boiler & Quartz Condenser & Quartz Heater with built in triple safety feature).



TECHNICAL SPECIFICATIONS:

CAT.NO:	EP-DQ-STILL 2	EP-DQ-STILL 4	EP-DQ-STILL 8
Output	2 litres/hr, Double distilled	4 litres/hr, Double distilled	8 litres/hr, Double distilled
Boiler	QUARTZ	QUARTZ	QUARTZ
Condenser	QUARTZ	QUARTZ	QUARTZ
pH	6.9-7.0	6.9-7.0	
Conductivity, μScm^{-1}	0.5-0.9	0.5-0.9	0.5-0.9
Temperature	50-60°C	50-60°C	50-60°C
Pyrogen content *	Pyrogen free	Pyrogen free	Pyrogen free
Water supply	2 litre/min x 2 Input	2 litre/min x 2 Input	4 litre/min x 2 Input
Electricity supply	220 or 240V, 50 -60Hz, single phase	220 or 240V, 50 -60Hz, double phase	220 or 240V, 50 -60Hz, three phase
Power requirement	3kW	6kW	12kW

The above output result can be achieved by using treated water as a boiler feed.