

QUARTZ DOUBLE DISTILLER (Demountable vertical unit) **EQDD PANEL SERIES (Vertical model).**



TECHNICAL SPECIFICATIONS:

TYPE	OUTPUT CAPACITY	POWER REQUIRED	COOLINWATERFEED (approx)
EQDD1	1ltrs/hr	2kw Single Phase	40Liters/hr
EQDD1.5	1.5 ltrs/hr	3kw Single Phase	60Liters/hr
EQDD2	2 ltrs/hr	4kw Single Phase	65Liters/hr
EQDD2.5	2.5 ltrs/hr	4.5kw Double Phase	70Liters/hr
EQDD3	3 ltrs/hr	5kw Double Phase	100Liters/hr
EQDD4	4 ltrs/hr	6kw Double Phase	120Liters/hr
EQDD5	5 ltrs/hr	8kw Double Phase	150Liters/hr

The above output result can be achieved by using treated water as a boiler feed.

POR* Price on Request