


QUARTZ DOUBLE DISTILLER EQDD-SHSERIES **(Demountable with Replaceable Silica Heater and Built in Triple Safety fail safe device)**



TECHNICAL SPECIFICATIONS:

TYPE	OUTPUT CAPACITY	POWER REQUIRED	COOLING WATER FEED
			(approx)
EQDD-SH1	1ltrs/hr	2kw	40Liters/hr
EQDD SH1.5	1.5 ltrs/hr	3kw	60Liters/hr
EQDD SH2	2 ltrs/hr	4kw	65Liters/hr
EQDD SH2.5	2.5 ltrs/hr	4.5kw	70Liters/hr
EQDD SH3	3 ltrs/hr	5kw	100Liters/hr
EQDD SH4	4 ltrs/hr	6kw	120Liters/hr
EQDD SH5	5 ltrs/hr	8kw	150Liters/hr

The above output result can be achieved by using treated water as a boiler feed.