


Single distillation with Quartz Boiler and Borosilicate Condenser **With built in triple safety features.**



TECHNICAL SPECIFICATIONS:

CAT.NO:	EAQ-2S-QB	EAQ-4S-QB	EAQ-8S-QB
Output	2 litres/hr, single distilled	4 litres/hr, single distilled	8 litres/hr, single distilled
Boiler	Quartz	Quartz	Quartz
Condenser	Borosilicate	Borosilicate	Borosilicate
Ph	5.5 - 6.0	5.5 - 6.0	5.5 - 6.0
Conductivity, μScm^{-1}	1.0-2.0	1.0-2.0	1.0-2.0
Temperature	30 - 40°C	30 - 40°C	30 - 40°C
Pyrogen content *	Pyrogen free	Pyrogen free	Pyrogen free
Water supply	2 litre/min	2 litre/min	4 litre/min
Electricity supply	220 or 240V, 50-60Hz, single phase	220 or 240V, 50-60Hz, single phase	220 or 240V, 50-60Hz, double phase
Power requirement	1.5kW	3kW	6kW

The above output result can be achieved by using treated water as a boiler feed.