


## **Quartz Single Distillation.** **With Built in Triple Safety Features**



### **TECHNICAL SPECIFICATIONS:**

<b>CAT.NO:</b>	<b>EAQ-2SQ</b>	<b>EAQ-4SQ</b>	<b>EAQ-8SQ</b>
Output	2 litres/hr, single distilled	4 litres/hr, single distilled	8 litres/hr, single distilled
Boiler	Quartz	Quartz	Quartz
Condenser	Quartz	Quartz	Quartz
pH	5.0-8.0	5.0-8.0	5.0-8.0
Conductivity, $\mu\text{Scm}^{-1}$	0.5-0.8	0.5-0.8	0.5-0.8
Temperature	30 - 40°C	30 - 40°C	30 - 40°C
Pyrogen content *	Pyrogen free	Pyrogen free	Pyrogen free
Water supply	2 litre/min	2 litre/min	4 litre/min
Electricity supply	220 or 240V, 50-60Hz, single phase	220 or 240V, 50-60Hz, single phase	220 or 240V, 50-60Hz, double phase
Power requirement	1.5kW	3kW	6kW

The above output result can be achieved by using treated water as a boiler feed.