

QUARTZ SINGLE DISTILLATION (Quartz Boiler & Quartz Condenser) EQSD PANEL SERIES (Vertical model).



TECHNICAL SPECIFICATIONS:

TYPE	OUTPUT CAPACITY	POWER REQUIRED	COOLING WATER FEED	BIOLOGICAL PURITY	RH/PH	CONDUCTIVITY	DISTILLATE TEMPERATURE
EQSDQ1	1.0	220-240V	30Liters/hr	PYROGEN FREE	6.0-7.0	1.5-2.0	60° ± 5° c
		1kwSingle Phase	(approx)				
EQSDQ2	2.0	220-240V	40Liters/hr (approx)	PYROGEN FREE	6.0-7.0	1.5-2.0	60° ± 5° c
		1.5kwSingle Phase					
EQSDQ2.5	2.5	220-240V	50Liters/hr (approx)	PYROGEN FREE	6.0-7.0	1.5-2.0	60° ± 5° c
		1.75kwSingle Phase					
EQSDQ3	3.0	220-240V	55Liters/hr	PYROGEN FREE	6.0-7.0	1.5-2.0	60° ± 5° c
		2kwSingle Phase	(approx)				
EQSDQ4	4.0	220-240V	60Liters/hr (approx)	PYROGEN FREE	6.0-7.0	1.5-2.0	60° ± 5° c
		3kwSingle Phase					
EQSDQ5	5.0	220-240V 4kw	65Liters/hr (approx)	PYROGEN FREE	6.0-7.0	1.5-2.0	60° ± 5° c
		Double Phase					
EQSDQ10	10.0	220-240V 8kw	120Liters/hr (approx)	PYROGEN FREE	6.0-7.0	1.5-2.0	60° ± 5° c
		Double Phase					

The above output result can be achieved by using treated water as a boiler feed.